



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse PHLPP (KIAA0606) Polyclonal Antibody, Rabbit

Mouse KIAA0606 (mKIAA0606) protein is a homologue of human KIAA0606 (ref. 1). KIAA0606 is identical to PH domain and leucine rich repeat protein phosphatase (PHLPP). Rabbit anti-mouse PHLPP (mPHLPP, mKIAA0606) antibody is raised against GST-fused recombinant protein (GX0356) containing following sequence:

(118 amino acids)

GSRVEVEVDIHCSRAKEKERQQHLLQVPAEASDEGIVISANEDESGLSKKADFSAVGTIGRRRANG
SVAPQE RSHNVIEVAADAPLRKPGGYFAAPAQPDQFIIPPELEEEVKEI

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0356. This antibody detects mPHLPP protein. It also recognizes human PHLPP protein. Details are shown in the InGaP database. (<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0606AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619