



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ACIN1 (KIAA0670) Polyclonal Antibody, Rabbit

Mouse KIAA0670 (mKIAA0670) protein is a homologue of human KIAA0670 (ref. 1). KIAA0670 is identical to Apoptotic chromatin condensation inducer in the nucleus (ACIN1). Rabbit anti-mouse ACIN1 (mACIN1, mKIAA0670) antibody is raised against GST-fused recombinant protein (GX0829) containing following sequence:

(105 amino acids)

FKRKISVVSATKGQAGNSDTEGGQPGRKRRWGASTAATQKKPSISITTESLKEAVVDLHADDSSR
ISEDETE RNGDDGTHDKGLKICRTVTQVVP AEGQENGQRE

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0829. This antibody detects mACIN1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA0670AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619