



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse TSPYL4 (KIAA0721) Polyclonal Antibody, Rabbit

Mouse KIAA0721 (mKIAA0721) protein is a homologue of human KIAA0721 (ref. 1). KIAA0721 is identical to TSPY-like 4 (TSPYL4). Rabbit anti-mouse TSPYL4 (mTSPYL4, mKIAA0721) antibody is raised against GST-fused recombinant protein (GX1772) containing following sequence:

(251 amino acids)

TGKEGEAGAAMQEKKGLQKEKKVAGGGKEETRPRAPKINCMSLEAIDQELSNVNAQADRAFLQL
ERKFGRMRRLHMQRSSFIIQNIPGFVWTAFRNHPQLSPMISGQDEDMRYMINLEVEELKQPRVG
CKFKFIFQSNPYFRNEGLVKEYERRSSGRVSVSLSTPIRWHRGQEPQAHHRNREGNTIPSFFNWFS
DHSLLFEFDRIAEIIGEL WSNPLQYYLMGDGPRRGVVRVPPRQPVESPRSFRRFQSG

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: mTSPYL4 protein. It also recognizes human TSPYL4 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0721AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619