



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse NRXN3 (KIAA0743) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0743PA

Gene: mouse neurexin 3 (mNRXN3, mKIAA0743)

Immunogen: GX1111 (GST-fusion protein, 353 amino acids)

VTFKTKSSYLTLATLQAYTSMHLFFQFKTTSADGFILFNSGDGNDFIAVELVKGYIHVFDLGNG
PNVIKGNDRPLNDNQWHNVITRDSSNTHSLKVDTKVVTQVINGAKNLDLKGDLYMAGLAQG
MYSNLPKLVASRDGFQGGCLASVDLNGRLPDINDALHRSGQIERGCEVALTKADLQGPSTTCQ
EDSCANQGVCMQQWEGFTCDCSMTSYSGNQCNPDGATYIFGKSGGLILYTWPANDRPSTRS
DRLAVGFSTTVKDGILVRIDSAPGLGDFLQLHIEQGKIGVVFNIGTVDISIKEERTPVNDGKYHV
RFRTRNGGNATLQVDNWPVNEHYPTGRQLTIFNTQA
ROUGE Database (Mouse KIAA/FLJ cDNA Database)
[mKIAA0743](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1111. This antibody detects mNRXN3 protein. It also recognizes human NRXN3 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619