



COSMO BIO Co., LTD.  
Inspiration for Life Science

## Anti-Mouse RUFY3 (KIAA0871) Polyclonal Antibody, Rabbit

**Product Code:** PRX-MKA0871

**Gene:** mouse RUN and FYVE domain containing 3 (mRUFY3, mKIAA0871)

**Immunogen:** GX0978 (GST-fusion protein, 191 amino acids)

NNLQTKVDLLEKSNTKLTEELAVANNRIITLQEEMERVKEESSYLLESNRKGPQDRTAEGQAL

SEARKHLKEETQLRLDVEKELELQISMQRQEMELAMKMLEKDVCEKQDALVSLRQQQLDDLRLALK

HELAFLQSSDLGVKQKSELNSRLEEKTNQMAATIKQLEQSEKDLVKQAKTLNSAANKLIPKHH

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

**Format:** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity:** 100 µL

**Specificity:** Specific to recombinant protein GX0978. This antibody detects mRUFY3 protein. It also recognizes human RUFY3 protein. Other species have not been tested.

**Application:** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage:** Store at -20°C. Avoid freeze-thaw cycles

**References** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619