



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse FAM169A (KIAA0888) Polyclonal Antibody, Rabbit

Product Code: MKA0888

Gene: mouse Family with sequence similarity 169, member A (mFAM169A, mKIAA0888)

Immunogen: GX2302 (GST-fusion protein, 147 amino acids)

KKTLTGSSDNVATVSNIEKSDGNFNSVVPFPEEPVSQNLSPNTTSSVEDQGEEGAPEAQEP
SATQSSLIEVELEDAPFPQAGQKSQSEEQSEASSEHLEQFTQSAEKAVDSSSEEIEVEVPVVD
RRNLRRKAKGHKGPGKKKAK

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA0888](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX2302. This antibody detects mFAM169A protein. It also recognizes human FAM169A protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619