



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse USP19 (KIAA0891) Polyclonal Antibody, Rabbit

Code # : MKA0891PA [FIG01614], (Lot# XXX)

Gene : mouse ubiquitin specific peptidase 19 (mUSP19, mKIAA0891)

Immunogen : GX0401 (GST-fusion protein, 237 amino acids)
GYAQHDAQEFMAFLDGLHEDLNRIQNKPYTETVDS DGRPDEVVAEEAWQRHKMRNDSFIVD
LFQGQYKSKLVCPVCAKVSITFDPLYLPVPLPQKQKVLPIFYFAREPHSKPIKFLVSVSKENSSA
SEVLDSLQSVHVKPENLRRLAEVIKNRFRVFLPSHSLDAVSPTDVLFCFELLSPELAKERVVVL
EVQQRQVPSIPISKCAACQQRKQQSEEEKLRCTRCYRVGYCNQ

Format : Rabbit IgG purified with Protein A affinity chromatography.

Constitution : PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity : 100 µg (200 µL)

Specificity : Specific to recombinant protein GX0401. This antibody detects mUSP19 protein. It also recognizes human USP19 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619