



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse TRIM37 (KIAA0898) Polyclonal Antibody, Rabbit

Mouse KIAA0898 (mKIAA0898) protein is a homologue of human KIAA0898 (ref. 1). KIAA0898 is identical to tripartite motif-containing protein 37 (Mulibrey nanism protein, TRIM37). Rabbit anti-mouse TRIM37 (mTRIM37, mKIAA0898) antibody is raised against GST-fused recombinant protein (GX0196) containing following sequence:

(254 amino acids)

STDMQTNLFCADQAALTTCGPENSGRLQDLGMELLAKSSVAGCYIRNPTNKKNSPKSARAIAGSLS
LRRVDSGENSRSKGDCQVLAEGSSGSSQSGSRHSSPRALTHGIIGDLLPKSEDRQCKALDSAV
VVAVFNGLPVTEKRRKMVTLGTNAKGGRLEGMMQADLESHSEAGEVQPTLPEGASAAPEEGMSS
DSDIECDTENEEQEEHT SMGAFNDPFLAQPPDEDSHSSFPDGEQIDPENLHFNPDDEGGG

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0196. This antibody detects mTRIM37 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Immunoprecipitation.

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0898AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619