



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse CAMTA2 (KIAA0909) Polyclonal Antibody, Rabbit

Mouse KIAA0909 (mKIAA0909) protein is a homologue of human KIAA0909 (ref. 1). KIAA0909 is identical to calmodulin binding transcription activator 2 (CAMTA2). Rabbit anti-mouse CAMTA2 (mCAMTA2, mKIAA0909) antibody is raised against GST-fused recombinant protein (GX0828) containing following sequence:

(202 amino acids)

AAQGYARLIETLSQWRSVETGSLDLEQEVDP LNVDHFSCTPLMWACALGHLEAAVLLFCWNRQALSI
PDSLGRPLSVAHSGRHVRLARCLEELQRQELSVEHPLALSPSSSPDTGLSSASSPSELSDGTFVS
TSAYSSAPDGSPPPPAPPLASDISMEMIPGQLSCGAPETPLLLMDYEATNSKEPAPSPCGPPLAQDNG
A

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0828. This antibody detects mCAMTA2 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0909AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619