



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse NAV3 (KIAA0938) Polyclonal Antibody, Rabbit

Mouse KIAA0938 (mKIAA0938) protein is a homologue of human KIAA0938 (ref. 1). KIAA0938 is identical to Neuron navigator 3 (Steerin-3, Pore membrane and/or filamentinteractinglike protein 1, Unc-53 homolog 3, NAV3). Rabbit anti-mouse NAV3 (mNAV3, mKIAA0938) antibody is raised against GST-fused recombinant protein (GX0145) containing following sequence:

(143 amino acids)

FLPCPMDVEGSRVWFMDLWNYSLVYPVLEAVREGLQMYGKRAPWEDPSKWVLDTPWSSASL
PQEGPALLQLRPEDVGYEACTSTKEATTSKHIPQTDTEGDPLMNMLMKLQEAANYPSTQSCDG
DSVSHREDILDTSES TL

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0145. This antibody detects mNAV3 protein. Other species have not been tested.

Application: Western blotting

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0938AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619