



COSMO BIO Co., LTD.  
Inspiration for Life Science

## Anti-Mouse TBC1D30 (KIAA0984) Polyclonal Antibody, Rabbit

**Product Code:** MKA0984

**Gene:** mouse fibronectin type III domain containing 3A (mTBC1D30, mKIAA0984)

**Immunogen:** GX1640 (GST-fusion protein, 278 amino acids)

FEGSEIILRVSLAIWAKLGEQIECCETAEDEFYGTMGRLTQEMLEQDLLQSHELMQTVYSMAFPF  
FPQLAELREKYTYNITPFPATIKPTSVSGRHSKARDSDDENGPDDEDASAVGCLGPLSGLLA  
PELQKYQKQIKEATEEQTLRSNNIAELSPGAINSCRSEYHAAFNSMMMERMTTDINALKRQYSRI  
KKKQQQLHQVYIRADKGPVTSILPSQANSSPVINHLLLGKKMKITNRAAKNAVIHVP GHPGGKI  
SPVPYEDIKTLNSPWRTHI

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA0984](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

**Format:** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity:** 100 µL

**Specificity:** Specific to recombinant protein GX1640. This antibody detects mTBC1D30 protein. It also recognizes human TBC1D30 protein. Other species have not been tested.

**Application:** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage:** Store at -20°C. Avoid freeze-thaw cycles

**References** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619