



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse CDC2L6 (KIAA1028) Polyclonal Antibody, Rabbit

Product Code: MKA1028

Gene: mouse Cell division cycle 2-like protein kinase 6, CDC2- related protein kinase 6, Death-preventing kinase, Cyclin-dependent kinase 11 (mCDC2L6, mKIAA1028)

Immunogen: GX0409 (GST-fusion protein, 302 amino acids)

EHDLWHIIKFHRASKANKKPMQLPRSMVKSLLYQILDGIHYLHANWVLHRDLKPANILVMGEGP
ERGRVKIADMGFARLFNSPLKPLADLPVVVTFWYRAPELLLGARHYTKAIDIWAIGCIFAELLTS
EPIFHCRQEDIKTSNPFHHDQLDRIFSVMGFPADKDWEDIRKMPEYPTLQKDFRRTTYANSSLIK
YMEKHVKVPDSKVFLLLQKLLTMDPTKRITSEQALQDPYFQEDPLPLTDVFAGCQIPYPKREFLN
EDEPEEKGDKNQPPQQNPHQQPAAPQAAPPQAPPQSS
ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA1028](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX0409. This antibody detects mCDC2L6 protein. It also recognizes human CDC2L6 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619