



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KCND2 (KIAA1044) Polyclonal Antibody, Rabbit

Mouse KIAA1044 (mKIAA1044) protein is a homologue of human KIAA1044 (ref. 1). KIAA1044 is identical to zinc finger RNA binding protein 2 (KCND2). Rabbit anti-mouse KCND2 (mKCND2, mKIAA1044) antibody is raised against GST-fused recombinant protein (GXD01) containing following sequence:

(187 amino acids)

YMQSKRNGLLSNQLQSSSEDEPAFISKSGSSFETQHHLLHCLEKTTNHEFVDEQVFEESCMEVAT
VNRPSHRPSLSSQQGVSTCCSRRHKKTFRIPNANVSGSHRGSVQELSTIQIRCVERTPLSNSRS
SLNAKMEECVKL NCEQPYVTTAIIPTPPVTTPEGDDRPEPEYSGGNIVRVSA

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GXD01. This antibody detects mKCND2 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1044AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619