



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ERC1 (KIAA1081) Polyclonal Antibody, Rabbit

Mouse KIAA1081 (mKIAA1081) protein is a homologue of human KIAA1081 (ref. 1). KIAA1081 is identical to ELKS/RAB6-interacting/CAST family member 1 (ERC1). Rabbit anti-mouse ERC1 (mERC1, mKIAA1081) antibody is raised against GST-fused recombinant protein (GX2134) containing following sequence:

(192 amino acids)

LTSRQVKDQNKVKANLKHKEQVEKKKSAQMLEEARREDLSLSDSSQQLQVEELLMAMEKVKQE
LESMKAKLSSTQQSLAEKETHLTNLRAERRKHLEEVLEMKQEALLAAISEKDANIALLELSSSKKK
TQEEVAALKREKDRVLVQQLKQQTQNRMKLMADNYEDDHFSSRSNQTNHKPSPDQDEEEGIW
A

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX2134. This antibody detects mERC1 protein. It also recognizes human ERC1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Immunoprecipitation.

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1081AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress

Manufactured by

ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619