



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse TRIM35 (KIAA1098) Polyclonal Antibody, Rabbit

Code # : MKA1098PA [FIG01625], (Lot# XXX)

Gene : mouse tripartite motif-containing 35 (mTRIM35, mKIAA1098)

Immunogen : GX1416 (GST-fusion protein, 149 amino acids)
AAPAMEPGPSVSPGPSRSFKEELLCAVCYDPFRDAVTLRCGHNFCRRCSVSGCWEVQTTPSCP
VCKERA VPGELRTNHTLNNLVETLLREEAEGARWTGRRSPRPCRAHRAPLTLFCLEDKELLCC
ACQADARHQEHRVQPIKDTAQDFR

Format : Rabbit IgG purified with Protein A affinity chromatography.

Constitution : PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity : 100 µg (200 µL)

Specificity : Specific to recombinant protein GX1416. This antibody detects mTRIM35 protein. It also recognizes human TRIM35 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619