



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse MYT1L (KIAA1106) Polyclonal Antibody, Rabbit

Mouse KIAA1106 (mKIAA1106) protein is a homologue of human KIAA1106 (ref. 1). KIAA1106 is identical to myelin transcription factor 1-like (MYT1L). Rabbit anti-mouse MYT1L (mMYT1L, mKIAA1106) antibody is raised against GST-fused recombinant protein (GX0194) containing following sequence:

(154 amino acids)

RATSAMKKAKLSGEQMLTIKQRASNGIENDEEIKQLDEEIKELNESNSQMEADMIKLRTQITTMESNLK
TIEEENKVIEQQNESLLHELANLSQSLIHSLANIQLPHMDPINEQNFDAYVTTLTEMYTNQDRYQSPEN
KALLENIKQAVRGIQV

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0194. This antibody detects mMYT1L protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA1106AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619