



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse REXO1 (KIAA1138) Polyclonal Antibody, Rabbit

Mouse KIAA1138 (mKIAA1138) protein is a homologue of human KIAA1138 (ref. 1). KIAA1138 is identical to RNA exonuclease 1 (REX1) homolog (REXO1). Rabbit anti-mouse REXO1 (mREXO1, mKIAA1138) antibody is raised against GST-fused recombinant protein (GX1676) containing following sequence:

(241 amino acids)

LVSSSGRCVRSEECYYHWGRLRRNRVAGGWETQYMCCSAVGSVGCQVAKQHVDGRKENLE
GFVRTFQKELPEDAHAGVFALDCEMSYYTYGLELTRVTVDMDMQVVDYDFVKPDNEVDYNTFR
SGVTEADLVDTSLTRDVQAVLLSMFSADTILIGHSLESDLLALKVIHGTVDTSVLFPHRLGLPYKRS
LRNLMADYLRQIIQDNVD GHSSSEDASACMHLVIWKIREDAKTKR

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1676. This antibody detects mREXO1 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1138AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619