



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ZNF512B (KIAA1196) Polyclonal Antibody, Rabbit

Mouse KIAA1196 (mKIAA1196) protein is a homologue of human KIAA1196 (ref. 1). KIAA1196 is identical to zinc finger protein 512B (ZNF512B). Rabbit anti-mouse ZNF512B (mZNF512B, mKIAA1196) antibody is raised against GST-fused recombinant protein (GX2373) containing following sequence:

(202 amino acids)

LLGVERTPSGRIRRTSAQVAVFHLQEIAEDELARDWTKRRMKDDLVPETARLNYTRPGLPTLNPQLLEA
WKNEVKEKGHVNCNPNDCCCEAIYASVSGLKAHLASCSKGDHLVGKYRCLLCPKEFSSESGVKYHILKTH
GENWFRTSADPSSKHKSQDSLMPRKEKKKSLSGGKKRGRKPKERSSEEPASKLPPNRDDWPPGGR

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX2373. This antibody detects mZNF512B protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1196AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619