



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse YEATS2 (KIAA1197) Polyclonal Antibody, Rabbit

Mouse KIAA1197 (mKIAA1197) protein is a homologue of human KIAA1197 (ref. 1). KIAA1197 is identical to YEATS domain containing 2 (YEATS2). Rabbit anti-mouse YEATS2 (mYEATS2, mKIAA1197) antibody is raised against GST-fused recombinant protein (GX0258) containing following sequence:

(211 amino acids)

GLKTFDPMFNFHPAIAKKFLESPSRSSSPTNQRSETPSANHSESDLSQHNDFLSDKDNNNSNVDVEE
RPPSTGEQRPSRKDTSSISGSHKRELRNADLTGDETSRLFVKKTIIVGVNSKYIPDPKREENDQSTH
KWMVYVRGSRREPSINHFVKKVWFFLHPSYKPNDLVEVREPPFHLTRRGWGEFPVRVQVHFKDSQ
NKRIDIIHNLKVL

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0258. This antibody detects mYEATS2 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA1197AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619