



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse GRAMD1B (KIAA1201) Polyclonal Antibody, Rabbit

Mouse KIAA1201 (mKIAA1201) protein is a homologue of human KIAA1201 (ref. 1). KIAA1201 is identical to GRAM domain containing 1B (GRAMD1B). Rabbit anti-mouse GRAMD1B (mGRAMD1B, mKIAA1201) antibody is raised against GST-fused recombinant protein (GX0439) containing following sequence:

(370 amino acids)

CEEIPIEENEVDSSSKSSIETKPDASPQLPKKSITNSTLTSTGSSEAPVSFDGLPLEEEVMEGDGSLEKEL
AIDNIIGEKIEIMAPVTSPSLDFNDNEDIPTELSDSSDTHDEGEVQAFYEDLSGRQYVNEVFNFSDKLYDL
LFTNSPFLRDFMEQRRFSDIIFHPWKKEENGNSRVILYTTITLNPLAPKTATVRETQTMKASQESECYVI
DAEVLTHDVPYHDYFYTINRYTLTRVARNKSRRLRVSTELRYRKQPWGFVKTFIEKNFWSGLEDYFRHLET
ELTKTESTYLAEIHRQSPKEKASKSSAVRRRRKRPHAHLRVPHLEEVMSPTTPTDEDVGHRIKHVAGSTQ
TRHIPEDTP DGFHL

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0439. This antibody detects mGRAMD1B protein. It also recognizes human GRAMD1B protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA1201AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619