



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse MYOF (KIAA1207) Polyclonal Antibody, Rabbit

Product Code: MKA1207PA

Gene: mouse Myoferlin (Fer-1-like protein 3) (mMYOF, mKIAA1207)

Immunogen: GX1489 (GST-fusion protein, 273 amino acids)

SLVVECGGERVESVVIKSLKKTNPFPSSVLFMKVFLPKEELYMPPLVIKVIDHRQFGRKPVVGGQC
TIDHLDRFRCDPYAGKEDIVPQLKASLMSAPPCREVVIEIEDTKPLLASKLSEKEEEIVDWWSKF
YASSGEHEKCGQYIQKGYSKLIKIDCELEDVADFEGLTDFSDTFKLYRGKSDENEDPSVVGFEFK
GSFRIYPLDDPSVPAPPRQFRELPSVPQECTVRIYIVQGLQLQPQDNNGLCDPYIKITLGKKVI
EDRDHYIPNTLNP

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA1207](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1489. This antibody detects mMYOF protein. It also recognizes human MYOF protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619