



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse CRSP3 (KIAA1216) Polyclonal Antibody, Rabbit

Product Code: MKA1216PA

Gene: mouse Cofactor required for Sp1 transcriptional activation subunit 3, Mediator of RNA polymerase II transcription subunit 23, Mediator complex subunit 23, CRSP complex subunit 3, Transcriptional coactivator CRSP130, Vitamin D3 receptor-interacting protein complex 130 kDa component, DRIP130, Activator-recruited cofactor 130 kDa component, ARC130, Sur-2 (mCRSP3, mKIAA1216)

Immunogen: GX1561 (GST-fusion protein, 257 amino acids)
DTYLKHAMNAREDNPWVPEDSYCKLIGRLVDTMAGKSPGPFPCDWRFEFPNPAHALHV
TCVELMALAVPGKDVGNALLNVVLKSQPLVPRENITAWMNAIGLIITLPEPYWIVLHDRVNVISS
SSLTSETEWVGYPFRLFDFTACHQSYSEMCSYTLALAHAVWHSSIGQLSLIPKFLTEALLPVV
KTEFQLLYVYHLVGPFLQRFQQRTRCMIEIGVAFYDMLLNVDQCSTHLNYMDPICDFLYHMK
ROUGE Database (Mouse KIAA/FLJ cDNA Database)
[mKIAA1216](#)
The mammalian expression vector of the full length human protein corresponding to the antibody is available!
[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1561. This antibody detects mCRSP3 protein. It also recognizes human CRSP3 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by  ProteinExpress Co., Ltd. and  Kazusa DNA Research Institute.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619