



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse NDRG2 (KIAA1248) Polyclonal Antibody, Rabbit

Mouse KIAA1248 (mKIAA1248) protein is a homologue of human KIAA1248 (ref. 1). KIAA1248 is identical to NDRG family member 2 (NDRG2). Rabbit anti-mouse NDRG2 (mNDRG2, mKIAA1248) antibody is raised against GST-fused recombinant protein (GX2322) containing following sequence:

(137 amino acids)

DPNAKGWMDWAAHKLTGLTSSIPDMILGHLFSQEELSGNSELIQKYRGIIQHAPNLENIELYWNSYN
NRRDLNFERGGETTLKCPVMLVVGDAQPHEDAVECNSKLDPTQTSFLKMADSGGQPQLTQPGK
LTEAFK

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX2322. This antibody detects mNDRG2 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1248AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619