



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse SH2B1 (KIAA1299) Polyclonal Antibody, Rabbit

Mouse KIAA1299 (mKIAA1299) protein is a homologue of human KIAA1299 (ref. 1). KIAA1299 is identical to SH2B adaptor protein 1 (SH2B1). Rabbit anti-mouse SH2B1 (mSH2B1, mKIAA1299) antibody is raised against GST-fused recombinant protein (GX2547) containing following sequence:

(171 amino acids)

ERWTHRFRERLRLSRGGGTLKDGAGMIQREELLSFMGAEEAAPDPAGVGRGGGAAGLTSGGGGQ
PQWQKCRLLLRSEGEggggsrleffvppkasrpRLSIPcSTITDVRTATALEMPDRENTFVVKVE
GPSEYILETSDAL HVKAWVSDIQECLSPGPCPAISPRPMTLPH

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX2547. This antibody detects mSH2B1 protein. It also recognizes human SH2B1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1299AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619