



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse NPLOC4 (KIAA1499) Polyclonal Antibody, Rabbit

Product Code: MKA1499PA

Gene: mouse Neuroligin-3 precursor (Gliotactin homolog) (mNLGN3, mKIAA1480)

Immunogen: GX1022 (GST-fusion protein, 205 amino acids)

FPSSLAGPSSEMETSTSVGLKAFGAPNVVEIDQYLSKQDGKIYRSRDPQLCRHGKPLGKCVH
CVPLEPFDEEDYLNHLEPPVKHMSFHAYIRKLTGGADKGKFVALENISCKIKSGCEGHLWPNGI
CTKCQPSAITLNRQKYRHVDNIMFENHTVADRFLDFWRKTGNQHFYLYGRYTEHKDIPLGIRA
EVAAYEPPQIGTQ

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA1499](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1022. This antibody detects mNPLOC4 protein. It also recognizes human NPLOC4 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619