



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse CHD8 (KIAA1564) Polyclonal Antibody, Rabbit

Mouse KIAA1564 (mKIAA1564) protein is a homologue of human KIAA1564 (ref. 1). KIAA1564 is identical to Chromodomain-helicase-DNA-binding protein 8 (CHD8). Rabbit anti-mouse CHD8 (mCHD8, mKIAA1564) antibody is raised against GST-fused recombinant protein (GX2228) containing following sequence:

(170 amino acids)

GFWFTELLERNQCFTSWVSNRPHCFWMTGFFNPQGFLTAMRQEITRANKGWALDNMVL CNE
VTKFMKDDISAPPTG VYVYGLYLEGAGWDKRNMKLIESKPKVLFELMPVIRIFAENNTARDPRL
YCCPIYKKPVRTDLN YIAAVDLKTAQAP EHWVLRGVALLCDVK

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX2228. This antibody detects mCHD8 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1564AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619