



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse MICALL1 (KIAA1668) Polyclonal Antibody, Rabbit

Product Code: MKA1668

Gene: mouse MICAL-like protein 1, Molecule interacting with Rab13, MIRab13 (mMICALL1, mKIAA1668)

Immunogen: GX1751 (GST-fusion protein, 292 amino acids)
LGQPSGEQMPQVRTKGSAGTHSTKPFSGATPTPFLLAGDRNPAPPVGSASPQLQIKSSCKENP
FNRKPSPSASPTVRKATKGAKPVRPPAPGHGFPLIKRKVQADQYIPEEDIYGEMDNIERQLDAL
EHSGVLLEEKLRGGANEGSEDDMLVDWFKLIHEKHLLVRRESELIYVFKQQNLEQRQADVEFEL
RCLLNKPEKDWTDREAREKVLQMQLMTLIEQRDAIVNCLDEDRQREEEEDKMLETMIKKKDF
QREAESDSKKKGKFKTIKVLKFLGNKREAKSKAPGDKS
ROUGE Database (Mouse KIAA/FLJ cDNA Database)
[mKIAA1668](#)
The mammalian expression vector of the full length human protein corresponding to the antibody is available!
[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX1751. This antibody detects mMICALL1 protein. It also recognizes human MICALL1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619