



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse MTMR12 (KIAA1682) Polyclonal Antibody, Rabbit

Product Code: MKA1682

Gene: mouse myotubularin related protein 12 (mMTMR12, mKIAA1682)

Immunogen: GX0851 (GST-fusion protein, 247 amino acids)

DTNMGRENLDQSKPLTLLTVWDWSVQFEPKAQTLLRNPLYVEKPKLDRGQQKGSRFKHQRQ
LSLPLTQSKSSPKRGFFRETDHLIKNLLGKRISKLINSSDDLQDSSREFYDSWHSKPTDYHGLL
LPHIEGPEIKVWAQRYLRWIPEAQILGGGRVATMGKLEMMEEVQSLQEKEARHHRQEAVHAE
APGLLRNSVRLSSLFPFALLQRHSAKPVLPSTSGWKALGGEDDLAKREDEFVDLGDV
ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA1682](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX0851. This antibody detects mMTMR12 protein. It also recognizes human MTMR12 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619