



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse LCOR (KIAA1795) Polyclonal Antibody, Rabbit

Mouse KIAA1795 (mKIAA1750) protein is a homologue of human KIAA1795 (ref. 1). KIAA1795 is identical to ligand dependent nuclear receptor corepressor (LCOR). Rabbit anti-mouse LCOR (mLCOR, mKIAA1795) antibody is raised against GST-fused recombinant protein (GX0472) containing following sequence:

(88 amino acids)

YRQYNSEILEEAISVVMMSGKMSVSKAQSIYGIPHSTLEYKVKERLGLTKNPPKKKMKLMRSEGPD
VSVKIELD PQGEAAQSANESKTE

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0472. This antibody detects mLCOR protein. It also recognizes human LCOR protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1795AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619