



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse PHF6 (KIAA1823) Polyclonal Antibody, Rabbit

Mouse KIAA1823 (mKIAA1823) protein is a homologue of human KIAA1823 (ref. 1). KIAA1823 is identical to PHD finger protein 6 (PHF6). Rabbit anti-mouse PHF6 (mPHF6, mKIAA1823) antibody is raised against GST-fused recombinant protein (GX0505) containing following sequence:

(81 amino acids)

SQPGATIGCEIKACVKTYHYHCGVQDKAKYIENMSRGIYKLYCKNHSGNDEEDEERESKSR
GRVAIDQ QLTQQQLNGN

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0505. This antibody detects mPHF6 protein. It also recognizes human PHF6 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1823AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619