



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse CCDC85A (KIAA1912) Polyclonal Antibody, Rabbit

Mouse KIAA1912 (mKIAA1912) protein is a homologue of human KIAA1912 (ref. 1). KIAA1912 is identical to coiled-coil domain containing 85A (CCDC85A). Rabbit anti-mouse CCDC85A (mCCDC85A, mKIAA1912) antibody is raised against GST-fused recombinant protein (GX1954) containing following sequence:

(157 amino acids)

VLSGSPEHFQKHRPGGSPEHTRHSGGSPEHLQKHALGGSLEHLPRARGTSPEHLKQHYGASPDHKKH
ASGSGSGGGSREGTLRRPAQEDSSSHHRNVYSGMNESTLSYVRQLEARVRQLEENRMLPQVWW
RKLGDAAGSCPGIRQHLSGNQYKGPM

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1954. This antibody detects mCCDC85A protein. It also recognizes human CCDC85A protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA1912AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619