



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse MEX3B (KIAA2009) Polyclonal Antibody, Rabbit

Product Code: MKA2009

Gene: mouse RNA-binding protein MEX3B, RING finger and KH domain-containing protein 3, RING finger protein 195 (mMEX3B, mKIAA2009)

Immunogen: GX1936 (GST-fusion protein, 190 amino acids)

PSRDKEPVFEVTGMPENVDRAREEIEAHIALRTGGIILTDENDFHANGTDVGFDLHHGSGGSG
PGSLWSKPTPSITPTPGRKPFSSYRNDSSSSLGSASTDSYFGGGTSGSAAATSLADYSPPSP
ALSFHNGNNNNNGNGYTYTAGEASVSPDGGPELQPTFDPAAPPPGTPLLWAQFERSPPGGG
ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA2009](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX1936. This antibody detects mMEX3B protein. It also recognizes human MEX3B protein. Other species have not been tested..

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619