

Anti- Mouse SMEK1 (KIAA2010) Polyclonal Antibody, Rabbit

Product Code: MKA2010

Gene: mouse SMEK homolog 1, suppressor of mek1 (Dictyostelium)
(mSMEK1, mKIAA2010)

Immunogen: GX1101 (GST-fusion protein, 172 amino acids)
ILRNHRYRRDARTLEDEEEMWFNTDEDDMEDGEAVVSPSDKTKNDDDIMDPISKFMERKKLKE
SEEKEVLLKTNLSGRQSPSKLSLSSGTKTNLTSQSSATSLPGSPGSPGSPGSPGSPGSPVPSK
TSQTAAITTKGGLVGLVDYPDDDEDDDEDEDKEDTLPVSKKAKFES
ROUGE Database (Mouse KIAA/FLJ cDNA Database)
[mKIAA2010](#)
The mammalian expression vector of the full length human protein corresponding
to the antibody is available!
[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 μ L

Specificity: Specific to recombinant protein GX1101. This antibody detects mSMEK1 protein. It also recognizes human SMEK1 protein. Other species have not been tested..

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619