



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ZNF354C (KIAA4218) Polyclonal Antibody, Rabbit

Mouse KIAA4218 (mKIAA4218) protein is a homologue of human zinc finger protein 354C (ZNF354C) (ref. 1). Rabbit anti-mouse ZNF354C (mZNF354C, mKIAA4218) antibody is raised against GST-fused recombinant protein (GX0838) containing following sequence:

(219 amino acids)

KPYKCSECEKAFGSSSTLIKHLRVHTGEKPYRCRQCQKAFSQCSTLTVHQRIHTGEKLYKCAECDKAF
NCRACLHRHQRIHTGEKPYKCAECGKGYSQFPSLAEHQRLHTGEQLCQCLQCGRFTFTRVSTLIEHQR
IHTGQKPYQCNECGKTFNQYSSFNEHRKIHTGEKLYTCEECKAFGCKSNLYRHQRIHTGEKPYQCN
QCGKAFSQYS FLTEHER

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0838. This antibody detects mZNF354C protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA4218AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619