



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse FAM13B (KIBB0031) Polyclonal Antibody, Rabbit

Product Code: MKB0031

Gene: mouse Family with sequence similarity 13, member B (mFAM13B, mKIBB0031)

Immunogen: GX2834 (GST-fusion protein, 168 amino acids)

ESGEAQLSPQAARMTHHPLGEDCPPVLSHRSLDFGQSQRFLHDPEALDFSSKALSFTIRRSS
FSSKDEKREDRTPYQLVKKLQKKIRQFEEQFERERNKPSYSIDIAANPKVLKWMTELTCLRKQI
KDAKHKNSDGEFAPQTRPRSNTLPKSFSSLDHEDGESEGE

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX2834. This antibody detects mFAM13B protein. It also recognizes human FAM13B protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619