



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse USP5 (KIBB3523) Polyclonal Antibody, Rabbit

Code # : MKB3523 [FSR01304], (Lot# XXX)

Gene : mouse ubiquitin specific peptidase 5 (isopeptidase T) (mUSP5, mKIBB3523)

Immunogen : GX2940 (GST-fusion protein, 164 amino acids)
AELSEEALLSVLPTIRVPKAGDRVHKDECAFSFDTPSEGGLYICMNTFLGFGKQYVERHFNKT
GQRVYLHLRRTTRRPKEEDTSAGTGDPPRKKPTRLAIGVEGGFDLTEDKFEFDEDVKIVILPDYLE
IARDGLGGLPDIVRDRVTSAVEALLSADSASRKQE

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2940. This antibody detects mUSP5 protein. It also recognizes human USP5 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619