



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human Prorenin Polyclonal Antibody, Rabbit

Product Code: PR40-PRN01SR

Gene: REN (renin precursor)

Immunogen: Synthetic peptide corresponding to the N-terminal part of prosegment of human prorenin.

Format: Rabbit Antiserum

Preservative 0.05 % NaN₃


Quantity: 100 µL

Specificity: This antibody specifically detects human prorenin. It do not crossreact with human renin. Other species have not been tested.

Application: Immunoblotting 1 : 2,000,
Other applications have not been tested.

Storage: - 20°C, Avoid freeze-thaw cycles.

PR40-PRN01SR



Western blot

P : Human Prorenin (10 ng)
R : Human Renin (10 ng)

For research use only, not for diagnostic use

Manufactured by



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619