



Anti-Human Salusin-beta Polyclonal Antibody, Rabbit

Product Code: SL01-SLB01AF

Gene: A splice variant of TORSIN 2A (TOR2A), Preprosalusin

Immunogen: Synthetic peptide corresponding to human Salusin-beta peptide.

Format: Affinity Purified Rabbit IgG

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 50 µg (100 µL)

Specificity: This antibody specifically detects human Salusin-beta. It does not crossreact with human preprosalusin and prosalusin. Other species have not been tested.

Application: Immunoblotting 1 : 2,000, ELISA, RIA.
Other applications have not been tested.

Storage: -20°C, Avoid freeze-thaw cycles.

References 1) Shichiri M. et al.: Nat Med, 9, 1166 (2003).
2) Sato K. et al.: Regul Pept., 158(1-3), 63 (2009).

For research use only, not for diagnostic use

Manufactured by



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619