

**Certificate of Analysis**

**Product:** DyLight™ 800 Conjugated Affinity Purified anti-Human IgG (H&L) [Goat] Minimum Cross Reactivity to Bovine, Chicken, Goat, Guinea Pig, Hamster, Horse, Mouse, Rabbit, Rat and Sheep Serum Proteins

**Code:** 609-145-123

**Lot #** 99999

**Size:** 100 µg

**Physical State:** Lyophilized

**Antibody Concentration:** 1.0 mg/ml (by UV absorbance at 280 nm)

**Label:** DyLight™ 800 (MW 1050)

**Fluorochrome/Protein Ratio:** 2.0 DyLight™ 800 per mole of Goat IgG

**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2


| Excitation Wavelength: | Emission Wavelength: |
|------------------------|----------------------|
| 770 nm (in PBS)        | 794 nm (in PBS)      |

**Stabilizer:** 10 mg/ml BSA, IgG and Protease free

**Preservative:** 0.01% (w/v) Sodium Azide

**Application Note(s):** This antibody is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This antibody is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

**Instrument compatibility:** The emission spectra for this DyLight™ conjugate match the principal output wavelengths of most common fluorescence instrumentation.



**Figure 1.** DyLight™ dyes can be used for two-color Western Blot detection with low background and high signal. Anti-tubulin was detected using a DyLight™ 680 conjugate. Anti-TNF $\alpha$  was detected using a DyLight™ 800 conjugate. The image was captured using the Odyssey® Infrared Imaging System developed by LI-COR.

|                                 |                                     |                |
|---------------------------------|-------------------------------------|----------------|
| <b>Recommended Dilution(s):</b> | <b>FLUORESCENT WESTERN BLOT</b>     | >1:20,000      |
|                                 | <b>FLISA</b>                        | >1:10,000      |
|                                 | <b>IMMUNOFLUORESCENT MICROSCOPY</b> | >1:5,000       |
|                                 | <b>OTHER APPLICATIONS</b>           | User Optimized |






**Storage Conditions:** Store vial at 4° C prior to reconstitution. Rehydrate with 0.1 ml of deionized water (or equivalent). This product is stable for several weeks at 4° C as an undiluted liquid. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Dilute only prior to immediate use. Expiration date is one (1) year from date of receipt if properly stored.

**Purity and Specificity:** This product was prepared from monospecific antiserum by immunoaffinity chromatography using Human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Human IgG and Human Serum. No reaction was observed against Bovine, Chicken, Guinea Pig, Hamster, Horse, Mouse, Rabbit, Rat and Sheep Serum Proteins. This antibody will react with heavy chains of Human IgG and with light chains of most Human immunoglobulins.

**Immunogen:** Human IgG, whole molecule

**Conjugation Chemistry:** N-hydroxysuccinimide (NHS) ester

## Properties of DyLight™ Fluorescent Dyes:


| Emission      | Color   | DyLight™ Dye | Ex/Em (nm)     | $\epsilon$ ( $M^{-1} cm^{-1}$ ) | Similar Dyes                       |
|---------------|---|--------------|----------------|---------------------------------|------------------------------------|
| Green         |  | 488          | 493/518        | 70,000                          | Alexa™ 488, Cy2®, FITC             |
| Yellow        |  | 549          | 550/568        | 150,000                         | Alexa™ 546, Alexa 555, Cy3®, TRITC |
| Red           |  | 649          | 646/674        | 250,000                         | Alexa™ 647, Cy5®                   |
| Near Infrared |  | 680          | 682/715        | 140,000                         | Alexa™ 680, Cy5.5®, IRDye™ 700     |
| Infrared      |  | <b>800</b>   | <b>770/794</b> | <b>270,000</b>                  | <b>IRDye™ 800</b>                  |

**Relevant Links:** See [www.rockland-inc.com/DyLight](http://www.rockland-inc.com/DyLight) for additional information. Assortments of DyLight™ conjugates can be ordered as 3 and 10 vial value packs.

## Related Products:

- #[611-141-122](#) DyLight™ 488 Conjugated Affinity Purified anti-RABBIT IgG (H&L) [GOAT] MX10
- #[611-142-122](#) DyLight™ 549 Conjugated Affinity Purified anti-RABBIT IgG (H&L) [GOAT] MX10
- #[611-143-122](#) DyLight™ 649 Conjugated Affinity Purified anti-RABBIT IgG (H&L) [GOAT] MX10
- #[611-144-122](#) DyLight™ 680 Conjugated Affinity Purified anti-RABBIT IgG (H&L) [GOAT] MX10
- #[611-145-122](#) DyLight™ 800 Conjugated Affinity Purified anti-RABBIT IgG (H&L) [GOAT] MX10
- #[MB-070](#) Blocking Buffer for Fluorescent Western Blotting

## DyLight™ 800 Fluorescence Spectra:



**Note:** This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. There is no implied license hereunder for any commercial use. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 326, Gilbertsville, Pennsylvania, USA.

DyLight™ is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. IRDye™ and Odyssey® are trademarks of LI-COR Inc. Cy2®, Cy3®, Cy5® and Cy5.5® are trademarks of GE Healthcare Inc. Alexa® is a trademark of Invitrogen Inc.