



Anti Human HGF Monoclonal Antibody reduced Alpha chain (Mouse)

Catalog No.:	SIM-2ZCKHH2-EX
Clone No.:	Hyb-T3121 The splenic lymphocytes from BALB/c mouse immunized with reduced ricombinant Human HGF were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T3121) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with reduced Human HGF alpha chain
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

