



Anti Human IL-2 Monoclonal Antibody (Mouse)

| | |
|---------------------|--|
| Catalog No.: | SIM-2ZCKIL2 |
| Clone No.: | Hyb-6512 The splenic lymphocytes from BALB/c mouse immunized with E. coli expression ricombinant Human IL-2 were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-6512) were screened. |
| Source: | BALB/C Mouse implantation ascites |
| Specificity: | This antibody specifically reacts with Human IL-2 |
| Purify: | Ammonium sulfate precipitation |
| Subclass: | IgG1 |
| Buffer: | Physiological salt solution containing 0.1% Sodium Azide |
| Quantity: | 0.5mL [0.5mg/mL] |
| Application: | ELISA |
| Storage: | Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles. |

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZCKIL2

抗ヒト インターロイキン 2 モノクローナル抗体(マウス)
Anti IL 2

コード : 2ZCKIL2

クローン番号 : Hyb- 6512

由来 : 大腸菌発現リコンビナント - ヒトインターロイキン (IL) 2 を
免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞
(NS-1) を融合して得たハイブリドーマを BALB/c マウスに
接種して得られた腹水

特異性 : ヒト IL-2 に特異的に反応

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL : (03) 5632-9610 FAX : (03) 5632-9619
e-mail : mail@cosmobio.co.jp
URL : <http://www.cosmobio.co.jp/>