



MONOCLONAL ANTIBODY

Catalog No. 2ZCKRH1

Anti Rat HGF Alpha chain

Product type	Primary antibodies
Host	Mouse
Source	BALB/c Mouse implantation ascites
Form	Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.5 mL
(Concentration)	(0.5 mg/mL)
Specificity	This antibody specifically reacts with Rat HGF alpha chain
Antigen	The splenic lymphocytes from BALB/c mouse immunized with recombinant Human HGF were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T6942) were screened.
Clone	Hyb-T6942
Isotype	IgG ₁

Application notes ELISA¹⁾, WB²⁾

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Rat, Human, Mouse (HGF alpha chain)

Storage Maintain at 2-8 degrees (-20 degrees for long term storage) in undiluted aliquots.

Avoid repeated freeze/thaw cycles.

References 1) Yamada A., et al: Biomedical Research 16: 105-114, 1995.

2) Hattori N, et al: American Journal of Pathology 164: 1091-1098, 2004

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619