



Anti Human Endothelin A receptor Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZETRA-EX
Clone No.:	Hyb-AH3A8 The splenic lymphocytes from BALB/c mouse immunized with E. coli expression Human Endothelin A C-terminus (Amino acid sequence 359-427) were fused to mouse myeloma cells (P3X63Ag8 U-1) and the hybrid cells (Hyb-AH3A8) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human Endothelin A
Purify:	Ammonium sulfate precipitation
Subclass:	IgG2a
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

