

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHB14

Anti HBs (Antigenic determinant w)

Product type	Primary antibodies
Host	Mouse
Source	Ascites
Form	Hybridoma cell line (Hyb-4111) was grown in ascitic fluid BALB/c Mouse Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	This antibody specifically reacts with HBs antigen (determinant w) The antibody reacts with subtypes (adw,ayw)
Antigen	Purified HBs antigen (Antigenic determinant w)
Clone	Hyb-4111 The splenic lymphocytes from BALB/c mouse immunized with purified HBs antigen were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-4111) were screened.
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.
Avoid repeated freeze/thaw cycles.

References

- 1) Okamoto H, et al: Hepatology 6: 354-359,1986.
- 2) Usuda S, et al: J.Immunol.Methods 87: 203-210, 1986.
- 3) Okamoto H, et al: J.Virology 61: 3030-3034, 1987.
- 4) Tachibana K, et al: Viral.Immunol. 2: 25-29, 1989.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHB14

抗 HBs (抗原決定基 w)

Product type	一次抗体
Host	マウス
Source	BALB/c マウス腹水 ハイブリドーマ (Hyb-4111) を BALB/c マウスに接種して得られた腹水
Form	硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	HBs 抗原 (抗原決定基 w) に特異的に反応 サブタイプ adw, ayw に反応
Antigen	精製 HBs 抗原
Clone	Hyb-4111 精製 HBs 抗原を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞(NS-1)を融合して得られたハイブリドーマ
Isotype	IgG ₁

Application notes

ELISA

Recommended use

Optimal dilutions/concentrations should be determined by the end user.

Recommended dilutions

Staining Pattern

Cross reactivity

Storage

2~8°C 3 ヶ月間
-20°C 3 年間

(凍結融解の繰り返しは避けてください)

References

- 1) Okamoto H, et al: Hepatology 6: 354-359, 1986.
- 2) Usuda S, et al: J.Immunol.Methods 87: 203-210, 1986.
- 3) Okamoto H, et al: J.Virology 61: 3030-3034, 1987.
- 4) Tachibana K, et al: Viral.Immunol. 2: 25-29, 1989.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619