



Anti Human Albumin Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHBAL2-EX
Clone No.:	Hyb-4761 The splenic lymphocytes from BALB/c mouse immunized With Human serum origin albumin were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-4761) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human Albumin (Monomer, Polymer) This antibody do not reacts with Bovine Albumin, Horse Albumin, Mouse Albumin and Rabbit Albumin
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

