



## **Anti Human Glycophorin B Monoclonal Antibody (Mouse)**

<b>Catalog No.:</b>	<b>SIM-2ZHBGPB-EX</b>
<b>Clone No.:</b>	<b>Hyb-5362C</b> <b>The splenic lymphocytes from BALB/c mouse immunized With Human erythrocyte were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-5632C) were screened.</b>
<b>Source:</b>	<b>BALB/C Mouse implantation ascites</b>
<b>Specificity:</b>	<b>This antibody specifically reacts with Human Glycophorin B (This antibody coheres with all of ABO blood group and MN blood group)</b>
<b>Purify:</b>	<b>Ammonium sulfate precipitation</b>
<b>Subclass:</b>	<b>IgM</b>
<b>Buffer:</b>	<b>Physiological salt solution containing 0.1% Sodium Azide</b>
<b>Quantity:</b>	<b>0.5mL [0.5mg/mL]</b>
<b>Storage:</b>	<b>Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.</b>

***For research use only; not for use as a diagnostic.***

