



Anti Human Hemoglobin Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHBHB1-EX
Clone No.:	Hyb-T2007 The splenic lymphocytes from BALB/c mouse immunized With Human erythrocyte origin hemoglobin were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T2007) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human Hemoglobin
Purify:	Ammonium sulfate precipitation and gel filtration
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

