



Anti HBc Monoclonal Antibody Antigenic determinant beta (Mouse)

Catalog No.:	SIM-2ZHC22-EX
Clone No.:	Hyb-3120 The splenic lymphocytes from BALB/c mouse immunized With purified HBc antigen were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-3120) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with HBc antigen (Antigenic determinant beta)
Purify:	Ammonium sulfate precipitation
Subclass:	IgG2a
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

