



## **Anti Human AFP Monoclonal Antibody Antigenic determinant a (Mouse)**

<b>Catalog No.:</b>	<b>SIM-2ZHCMA1-EX</b>
<b>Clone No.:</b>	<b>Hyb-2051</b> <b>The splenic lymphocytes from BALB/c mouse immunized With purified human alpha fetoprotein(AFP) were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-2051) were screened.</b>
<b>Source:</b>	<b>BALB/C Mouse implantation ascites</b>
<b>Specificity:</b>	<b>This antibody specifically reacts with Human AFP (Antigenic determinant a)</b>
<b>Purify:</b>	<b>Ammonium sulfate precipitation</b>
<b>Subclass:</b>	<b>IgG1</b>
<b>Buffer:</b>	<b>Physiological salt solution containing 0.1% Sodium Azide</b>
<b>Quantity:</b>	<b>0.5mL [0.5mg/mL]</b>
<b>Storage:</b>	<b>Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.</b>

***For research use only; not for use as a diagnostic.***

