



Anti Human HSP90 α β Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHSP20
Clone No.:	Hyb-K3720A The splenic lymphocytes from BALB/c mouse immunized with E.coli ricombinant human Heat-Shock Protein (HSP) 90 β were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-K3720A) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with common region of Human HSP 90α and Human HSP 90β (Amino acid residues 291-304)
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.4mL [1mg/mL]
Application:	Western Blot, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHSP20

抗ヒト HSP90 モノクローナル抗体(マウス)
Anti HSP90α

コード : 2ZHSP20

クローン番号 : Hyb- K3720A

由来 : 大腸菌発現リコンビナント - ヒト Heat-Shock Protein (HSP) 90 を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを BALB/c マウスに接種して得られた腹水

特異性 : ヒト HSP90 および の共通部位 (アミノ酸配列番号 ; 291-304) に特異的に反応

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 1.0 mg/mL

容量 : 0.4 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: <http://www.cosmobio.co.jp/>